



MYTHIC™ 22 AL

22 parameters haematology
analyzer Auto Loading system



SWISS HAEMATOLOGY
SOLUTION

Technical specifications

Sample volume	< 20 µl
Technology	<ul style="list-style-type: none"> • Impedance (WBC, RBC, PLT) • Spectrophotometry (HGB) • Optical (5-DIFF): Optical Cytometry Hydrofocusing Free Optical laser method for 5-part WBC differential using semiconductive diode
22 parameters	WBC, LYM% and #, MON% and #, NEU% and #, BAS% and #, EOS% and #, RBC, HGB, HCT, MCV, MCH, MCHC, RDW, PLT MPV, PCT, PDW
Curves	RBC curve, PLT curve
Scattergram	5-Diff scattergram
Dimensions	480 x 400 x 400 mm (W x D x H)
Weight	25 kg
Power	<ul style="list-style-type: none"> • 90 to 250 V • 50 to 60 Hz • Less than 150 VA
External printers	<ul style="list-style-type: none"> • Dot matrix printer • Monochrome inkjet printer, PCL3 type • Color laser printer, PCL6 type
User's interface	<ul style="list-style-type: none"> • TFT colour touchscreen 8"5 • Built-in numerical keyboard
Output	<ul style="list-style-type: none"> • Bidirectional LIS communication • RS232 • Ethernet (UDP/IP, TCP/IP protocol) • 5 USB ports
Tubes Id.	Internal barcode scanner
Loading	<ul style="list-style-type: none"> • Continuous loading • 10 racks of 5 tubes • STAT position • Automatic rerun option

Software specifications

Multilanguage	Bulgarian, Brazilian Portuguese, Chinese, English, French, Italian, Latvian, Polish, Romanian, Russian, Spanish						
Calibration	Automatic or manual						
Quality control	<ul style="list-style-type: none"> • 12 batches: 100 runs/batch • Table and Levey Jennings graphs 						
Archives	<ul style="list-style-type: none"> • Internal memory: 1500 patients with scattergram and histograms • Mass storage (USB): up to 60 000 patients 						
Traceability	<ul style="list-style-type: none"> • Patient identification: name, patient ID, sample ID, sex, ... • Daily log • Reagents monitoring: lot number, expiry date, volume • Service log 						
Flagging	<ul style="list-style-type: none"> • Pathologic flags • Lab limits • Reagents alerts • Instrument alerts 						
Reagents	<table border="0"> <tr> <td>HM 22-001-1</td> <td>• Cleaner for Mythic™ 22</td> </tr> <tr> <td>HM 22-002-1</td> <td>• OnlyOne™</td> </tr> <tr> <td>HM 22-003-10</td> <td>• Diluent for Mythic™ 22</td> </tr> </table>	HM 22-001-1	• Cleaner for Mythic™ 22	HM 22-002-1	• OnlyOne™	HM 22-003-10	• Diluent for Mythic™ 22
HM 22-001-1	• Cleaner for Mythic™ 22						
HM 22-002-1	• OnlyOne™						
HM 22-003-10	• Diluent for Mythic™ 22						
Calibrator	Myt-Cal: calibrator for Mythic™ 22						
Controls	Myt-5D tri level or normal level: quality controls for Mythic™ 22						



22 parameters haematology analyzer Auto Loading system

- 1 Orphée reserves own rights to technical and software modification. In case of printer choice please contact your local Orphée's representative.